



PGL-VI-561/1~5mW


**GREEN LASER MODULE
AT 561 nm**








PGL-VI green laser module at 561nm is made features of small size, high output power, long lifetime, which made it widely used in collimation, laser medical treatment, scientific experiment, optical instrument and so on.



SPECIFICATIONS

Wavelength (nm)	561±1
Operating mode	CW
Output power (mW)	>5
Transverse mode	TEM ₀₀
Beam diameter at the aperture (1/e ² , mm)	~1.5
Beam divergence, full angle (mrad)	<1.5
Operating temperature (°C)	15~30
Power supply	DC 3V PCB
TTL modulation	Optional, 1Hz-1KHz, 1KHz-10KHz, or higher on request
Expected lifetime (hours)	10000
Warranty	6 months
Note: The shell is positive electrodes.	



<p>PGL-VI-561</p>  <p>Φ20×100(L) mm², 0.1 kg</p>	<p>PGL-VI-1-561</p>   <p>Φ15×65(L) mm², 0.036 kg</p>	
 	<p>3VPCB</p>  <p>38×14×8(H) mm³</p>  <p>18×14×8(H) mm³</p>	